

# Professionelles Großhandels - Tourniquet für Erste -Hilfe- und Notfallversorgung für die medizinische Versorgung

## Produktdetail

|                    |                                                                                                                           |
|--------------------|---------------------------------------------------------------------------------------------------------------------------|
| <b>Produktname</b> | <b>Tourniquet</b>                                                                                                         |
| <b>Farbe</b>       | Schwarz                                                                                                                   |
| <b>Größe</b>       | 3,8*65 cm, 3,8*75 cm, 3,8*95 cm                                                                                           |
| <b>Material</b>    | Hakenschleifenband und Plastik/mentale Stange                                                                             |
| <b>MOQ</b>         | Kleine Bestellung akzeptieren                                                                                             |
| <b>Paket</b>       | 1 PC/OPP -Tasche oder individuelles Paket                                                                                 |
| <b>Anpassung</b>   | Das benutzerdefinierte Logo, die Verpackung und die Farbe können gemäß den Benutzeranforderungen erfolgen                 |
| <b>Funktion</b>    | Mit einer Hand angewendet, schnell und praktisch, zufällige Kratzer, Schnitte und andere gute Auswirkungen der Hämostase. |

**Tourniquet ist eine häufige lebensrettende Ausrüstung, die aus hochwertigen Materialien besteht und mit einer Hand betrieben werden kann, wodurch der Blutfluss vollständig an den Gliedmaßen durch schwere Blutungsverletzungen blockiert wird. Es ist ein großartiges Erste -Hilfe -Tool mit einer Vielzahl von Anwendungsszenarien, nicht nur in Kampf- und medizinischer Rettung, sondern auch in Außen-, Abenteuer- und arteriellen Blutungen. Ein Muss für alle!**

---

# CAT ROTARY TOURNIQUET



Easy to use



Breathable  
&  
Comfortable



Effective  
Hemostasis



Compact  
&  
Convenient





# FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



# HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response



# RECORD USAGE TIME

Write down the time of use to easily determine the length of time



# TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control



## PRODUCT USAGE



**Wilderness  
Combat**



**Outdoor  
Adventure**

## PREPAREDNESS AVERTS PERIL



**Bleeding  
from injury**



**Hiking and  
camping**

# — HOW TO USE —



**01**  
Tie the tourniquet three fingers' width above the wound to stop bleeding



**02**  
Apply tourniquet and tighten to stop bleeding



**03**  
Pull the hook loop tape tight and paste it back on



**04**  
Rotate the pressure fixed rod to tighten the strap



**05**  
Fasten enough to use the clasp to fix the rod



**06**  
Check for bleeding and distal pulse loss



**07**  
Attach the white hook loop tape



**08**  
Write down the time to prevent muscle necrosis caused by too much time

---

## Unternehmensinformationen



# CONTACT INFORMATION

TEL : +80 20 34992808  
+86 18026363923

Web : [www.homecommoditysupplier.com](http://www.homecommoditysupplier.com)

E-mail : [sales02@balance-china.com](mailto:sales02@balance-china.com)

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,  
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470