CAT rotierendes Tourniquet

Produktdetail

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Adjustable to fit arm and leg size

HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response

RECORD USAGE TIME

Write down the time of use to easily determine the length of time

TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control

PRODUCT USAGE

Wilderness
Combat

Outdoor
Adventure

PREPAREDNESS AVERTS PERIL

Bleeding
from injury

Hiking and
camping

— HOW TO USE —

01
Tie the tourniquet three fingers' width above the wound to stop bleeding

02
Apply tourniquet and tighten to stop bleeding

03
Pull the hook loop tape tight and paste it back on

04
Rotate the pressure fixed rod to tighten the strap

05
Fasten enough to use the clasp to fix the rod

06
Check for bleeding and distal pulse loss

07
Attach the white hook loop tape

08
Write down the time to prevent muscle necrosis caused by too much time

Firmeninformation

We have 5000m² factory, our own factory
The REAL STRENGTH is to show the true face to you

Strict certification testing provides quality guarantee and **WORRY-FREE EXPORT**

We have achieved some inspection certification, such as ISO9001:2015, BSCI, SGS, TUV, etc. and have won the title of “Guangdong Province Credible Enterprise” for three consecutive years. In addition, we have obtained over 20 product patents and trademark, as well as 5 patent evaluation reports. A series of quality inspection reports provide quality guarantee and worry-free export.

CONTACT INFORMATION

TEL : +80 20 34992808
+86 18026363923

Web : www.homecommoditysupplier.com

E-mail : sales02@balance-china.com

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470