

Il laccio d'ingrosso professionale per forniture mediche di pronto soccorso e di emergenza

Dettaglio del prodotto

Nome prodotto	Laccio
Colore	Nero
Misurare	3.8*65 cm, 3.8*75cm, 3.8*95cm
Materiale	Nastro anello di gancio e canna di plastica/mentale
Moq	Accetta un piccolo ordine
Pacchetto	1 PC/OPP Bag o Pacchetto personalizzato
Personalizzazione	Il logo, l'imballaggio e il colore personalizzati possono essere realizzati in base ai requisiti dell'utente
Funzione	Applicato con una mano, rapidi e pratici, graffi accidentali, tagli e altri buoni effetti dell'emostasi.

Il laccio emostatico è un'attrezzatura comune salvavita, che è realizzata con materiali di alta qualità e può essere gestito con una mano, bloccando completamente il flusso sanguigno agli arti da gravi ferite da trauma sanguinante. È un ottimo strumento di pronto soccorso con una vasta gamma di scenari di applicazione, non solo in combattimento e salvataggio medico, ma anche in sanguinamento all'aperto, avventura e arterioso. Un must per tutti!



CAT ROTARY TOURNIQUET



Easy to use



Breathable
&
Comfortable



Effective
Hemostasis



Compact
&
Convenient





Hook loop tape

TIME:

Time recording tape

Pressure fixed rod

Single Route Buckle

FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response



RECORD USAGE TIME

Write down the time of use to easily determine the length of time



TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control



PRODUCT USAGE



Wilderness
Combat



Outdoor
Adventure

PREPAREDNESS AVERTS PERIL



Bleeding
from injury



Hiking and
camping

— HOW TO USE —



01
Tie the tourniquet three fingers' width above the wound to stop bleeding



02
Apply tourniquet and tighten to stop bleeding



03
Pull the hook loop tape tight and paste it back on



04
Rotate the pressure fixed rod to tighten the strap



05
Fasten enough to use the clasp to fix the rod



06
Check for bleeding and distal pulse loss



07
Attach the white hook loop tape



08
Write down the time to prevent muscle necrosis caused by too much time

Informazioni sulla società

CONTACT INFORMATION

TEL : +80 20 34992808
+86 18026363923

Web : www.homecommoditysupplier.com

E-mail : sales02@balance-china.com

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470