

# Medical Tourniquet

## Product Detail

|                      |   |
|----------------------|---|
| <b>Product name</b>  | <b>CAT Rotary Tourniquet</b>  |
| <b>Color</b>         | <b>Black</b>  |
| <b>Size</b>          | <b>3.8*65cm, 3.8*75cm, 3.8*95cm</b>   |
| <b>MOQ</b>           | <b>Accept small order</b>   |
| <b>Package</b>       | <b>1 pc/opp bag,or customized package</b>   |
| <b>Customization</b> | <b>Customized logo, packaging and color can be made according to user requirements</b>                          |
| <b>Feature</b>       | <b>Applied with one hand, fast and practical, accidental scrapes, cuts and other good effect of hemostasis.</b> |

Tourniquet is a common life-saving equipment, which is made of high-quality materials and can be operated with one hand, completely blocking the blood flow to the limbs from severe bleeding trauma wounds. It is a great first aid tool with a wide range of application scenarios, not only in combat and medical rescue, but also in outdoor, adventure and arterial bleeding. A must-have for everyone!

# CAT ROTARY TOURNIQUET



Easy to use



Breathable  
&  
Comfortable



Effective  
Hemostasis



Compact  
&  
Convenient







# FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



## HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response





# RECORD USAGE TIME

Write down the time of use to easily determine the length of time



# TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control





## PRODUCT USAGE



**Wilderness  
Combat**



**Outdoor  
Adventure**

## PREPAREDNESS AVERTS PERIL



**Bleeding  
from injury**



**Hiking and  
camping**



# — HOW TO USE —



**01**  
Tie the tourniquet three fingers' width above the wound to stop bleeding



**02**  
Apply tourniquet and tighten to stop bleeding



**03**  
Pull the hook loop tape tight and paste it back on



**04**  
Rotate the pressure fixed rod to tighten the strap



**05**  
Fasten enough to use the clasp to fix the rod



**06**  
Check for bleeding and distal pulse loss



**07**  
Attach the white hook loop tape



**08**  
Write down the time to prevent muscle necrosis caused by too much time

## Company Information



**We have 5000m<sup>2</sup> factory, our own factory**  
**The REAL STRENGTH is to show the true face to you**





## Strict certification testing provides quality guarantee and **WORRY-FREE EXPORT**



We have achieved some inspection certification, such as ISO9001:2015, BSCI, SGS, TUV, etc. and have won the title of "Guangdong Province Credible Enterprise" for three consecutive years. In addition, we have obtained over 20 product patents and trademark, as well as 5 patent evaluation reports. A series of quality inspection reports provide quality guarantee and worry-free export.





## CONTACT INFORMATION

TEL : +80 20 34992808  
+86 18026363923

Web : [www.homecommoditysupplier.com](http://www.homecommoditysupplier.com)

E-mail : [sales02@balance-china.com](mailto:sales02@balance-china.com)

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,  
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470