

# Top Medical Tourniquet Supplier

## OEM Rotational Pressure

### Product Detail

Product name	Tourniquet
Color	Black
Size	3.8*65cm, 3.8*75cm, 3.8*95cm
Material	Hook loop tape & plastic/mental rod
MOQ	Accept small order
Package	1 pc/opp bag,or customized package
Customization	Customized logo, packaging and color can be made according to user requirements
Function	Applied with one hand, fast and practical, accidental scrapes, cuts and other good effect of hemostasis.

**Tourniquet is a common life-saving equipment, which is made of high-quality materials and can be operated with one hand, completely blocking the blood flow to the limbs from severe bleeding trauma wounds. It is a great first aid tool with a wide range of application scenarios, not only in combat and medical rescue, but also in outdoor, adventure and arterial bleeding. A must-have for everyone!**

---

# CAT ROTARY TOURNIQUET



Easy to use



Breathable  
&  
Comfortable



Effective  
Hemostasis



Compact  
&  
Convenient







# FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



# HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response





# RECORD USAGE TIME

Write down the time of use to easily determine the length of time



# TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control





## PRODUCT USAGE



**Wilderness  
Combat**



**Outdoor  
Adventure**

## PREPAREDNESS AVERTS PERIL



**Bleeding  
from injury**



**Hiking and  
camping**



# — HOW TO USE —



**01**  
Tie the tourniquet three fingers' width above the wound to stop bleeding



**02**  
Apply tourniquet and tighten to stop bleeding



**03**  
Pull the hook loop tape tight and paste it back on



**04**  
Rotate the pressure fixed rod to tighten the strap



**05**  
Fasten enough to use the clasp to fix the rod



**06**  
Check for bleeding and distal pulse loss



**07**  
Attach the white hook loop tape



**08**  
Write down the time to prevent muscle necrosis caused by too much time

## Company Information

### CONTACT INFORMATION

TEL : +86 20 34992808  
+86 18026363923

Web : [www.homecommoditysupplier.com](http://www.homecommoditysupplier.com)

E-mail : [sales02@balance-china.com](mailto:sales02@balance-china.com)

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,  
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470