

Torniquete por atacado profissional para fornecimento de suprimentos médicos de primeiros socorros e emergência

Detalhes do produto

Nome do produto	Tourniquet
Cor	Preto
Tamanho	3,8*65cm, 3,8*75cm, 3,8*95cm
Material	Fita de loop de gancho e haste de plástico/mental
MOQ	Aceitar um pequeno pedido
Pacote	1 PC/Opp Bag, ou pacote personalizado
Personalização	Logotipo personalizado, embalagem e cor podem ser feitos de acordo com os requisitos do usuário
Função	Aplicado com uma mão, rápida e prática, arranhões acidentais, cortes e outro bom efeito da hemostasia.

O torniquete é um equipamento comum que salva vidas, que é feito de materiais de alta qualidade e pode ser operado com uma mão, bloqueando completamente o fluxo sanguíneo para os membros de feridas de trauma de sangramento graves. É uma ótima ferramenta de primeiros socorros com uma ampla gama de cenários de aplicação, não apenas em combate e resgate médico, mas também em ao ar livre, aventura e sangramento arterial. Um must-have para todos!

CAT ROTARY TOURNIQUET



Easy to use



Breathable
&
Comfortable



Effective
Hemostasis



Compact
&
Convenient





FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response



RECORD USAGE TIME

Write down the time of use to easily determine the length of time



TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control



PRODUCT USAGE



**Wilderness
Combat**



**Outdoor
Adventure**

PREPAREDNESS AVERTS PERIL



**Bleeding
from injury**



**Hiking and
camping**

— HOW TO USE —



01
Tie the tourniquet three fingers' width above the wound to stop bleeding



02
Apply tourniquet and tighten to stop bleeding



03
Pull the hook loop tape tight and paste it back on



04
Rotate the pressure fixed rod to tighten the strap



05
Fasten enough to use the clasp to fix the rod



06
Check for bleeding and distal pulse loss



07
Attach the white hook loop tape



08
Write down the time to prevent muscle necrosis caused by too much time

Informações da empresa



CONTACT INFORMATION

TEL : +80 20 34992808
+86 18026363923

Web : www.homecommoditysupplier.com

E-mail : sales02@balance-china.com

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470