

Torniquete al por mayor profesional para primeros auxilios y suministros médicos de emergencia

Detalle del producto

Nombre del producto	Torniquete
Color	Negro
Tamaño	3.8*65cm, 3.8*75cm, 3.8*95cm
Material	Cinta de bucle de gancho y barra de plástico/mental
Moq	Acepta un pedido pequeño
Paquete	1 PC/Opp Bag, o paquete personalizado
Personalización	El logotipo personalizado, el embalaje y el color se pueden hacer de acuerdo con los requisitos del usuario
Función	Aplicado con una mano, raspados accidentales rápidos y prácticos, cortes y otros buenos efectos de la hemostasia.

El torniquete es un equipo común que salva la vida, que está hecho de materiales de alta calidad y se puede operar con una mano, bloqueando completamente el flujo sanguíneo a las extremidades de las heridas graves de trauma hemorragias. Es una gran herramienta de primeros auxilios con una amplia gama de escenarios de aplicación, no solo en combate y rescate médico, sino también en sangrado al aire libre, aventura y arterial. ¡Un imprescindible para todos!

CAT ROTARY TOURNIQUET



Easy to use



Breathable
&
Comfortable



Effective
Hemostasis



Compact
&
Convenient





FREE ADJUSTING BUCKLE

Adjustable to fit arm and leg size



HIGH-VISIBILITY RED TIP

In case of emergency, the red tip can immediately located and identified to make quick first aid response



RECORD USAGE TIME

Write down the time of use to easily determine the length of time



TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control



PRODUCT USAGE



**Wilderness
Combat**



**Outdoor
Adventure**

PREPAREDNESS AVERTS PERIL



**Bleeding
from injury**



**Hiking and
camping**

— HOW TO USE —



01
Tie the tourniquet three fingers' width above the wound to stop bleeding



02
Apply tourniquet and tighten to stop bleeding



03
Pull the hook loop tape tight and paste it back on



04
Rotate the pressure fixed rod to tighten the strap



05
Fasten enough to use the clasp to fix the rod



06
Check for bleeding and distal pulse loss



07
Attach the white hook loop tape



08
Write down the time to prevent muscle necrosis caused by too much time

Información de la empresa



CONTACT INFORMATION

TEL : +80 20 34992808
+86 18026363923

Web : www.homecommoditysupplier.com

E-mail : sales02@balance-china.com

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470