# Certified Tourniquet Wholesale Distributor Industrial Strength Tactical Manufacturer

#### **Product Detail**

| Product name  | Tourniquet   |
|---------------|--|
| Color         | Black  |
| Size          | 3.8*65cm, 3.8*75cm, 3.8*95cm   |
| Material      | Hook loop tape & plastic/mental rod  |
| MOQ           | Accept small order   |
| Package       | 1pc/opp bag,or customized package  |
| Customization | Customized logo, packaging and color can be made according to user requirements                          |
| Function      | Applied with one hand, fast and practical, accidental scrapes, cuts and other good effect of hemostasis. |

Tourniquet is a common life-saving equipment, which is made of high-quality materials and can be operated with one hand, completely blocking the blood flow to the limbs from severe bleeding trauma wounds. It is a great first aid tool with a wide range of application scenarios, not only in combat and medical rescue, but also in outdoor, adventure and arterial bleeding. A must-have for everyone!

## CAT ROTARY TOURNIQUET





## FREE ADJUSTING BUCKLE



### **HIGH-VISIBILITY RED TIP**

In case of emergency, the red tip can immediately located and identified to make quick first aid response



### **RECORD USAGE TIME**

Write down the time of use to easily determine the length of time



# TIGHTEN THE PRESSURE FIXED ROD

Rotatable internal band for quick bleeding control





#### **HOW TO USE-**



Tie the tourniquet three fingers' width above the wound to stop bleeding



Apply tourniquet and tighten to stop bleeding



Pull the hook loop tape tight and paste it back on



Rotate the pressure fixed rod to tighten the strap



Fasten enough to use the clasp to fix the rod



Check for bleeding and distal pulse loss



Attach the white hook loop tape



Write down the time to prevent muscle necrosis caused by too much time

#### **Company Information**

