

Dispositivo di selezione del catetere Adesivo impermeabile all'ordine di massa sterile Distributore di fornitori di prezzi di fabbrica

Deduzioni del prodotto

Ci sono innumerevoli usi per il gancio e il ciclo del porta del catetere con adesivo posteriore per le cure mediche. Il gancio e gli anelli afferrano bene in superficie e l'adesivo sul retro si tiene su qualsiasi superficie pulita e piatta. Tessuto di grado A, B e C, lodamento dell'adesione, un adesivo forte, l'adesione al freddo/resistente al calore e il diverso ciclo di gancio per la tua scelta! Siamo un produttore, forniamo prodotti di qualità e la migliore soluzione per i cusomeri "Progetto

Si prega di contattarci per maggiori dettagli.



- 1.Non-Woven Fabric
- 2.Strong Adhesive
- 3.Loop Side
- 4.Hook Side



Shape, size and color can be customized


※Please contact customer service for specific customization specifications

Widely used in medical field



- 1.Non-Woven Fabric**
- 2.Strong Adhesive**
- 3.Loop Side**
- 4.Hook Side**






※Contact customer service for more details



**15 YEARS
Of Manufacturing Experience**
Professional Manufacture Of Hook and Loop





Informazioni sulla società



CONTACT INFORMATION

TEL : +80 20 34992808
+86 18026363923

Web : www.homecommoditysupplier.com

E-mail : sales02@balance-china.com

Address : No.2 Shengyuan Street, Beiliu Ind Park, Dagang Town,
Nansha District, Guangzhou, Guangdong, China

Zip code: 511470